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Specimens required for Testing:

Hair samples:

In order to extract DNA from hair samples it is crucial to have hair roots. The roots are the white “swollen” end of the hairshaft, which is visible when plucked. The easiest and most reliable hair samples are found on the tail of cattle.

The following procedure should be followed:

1. The animal should be confined and easily able to be handled. Take hold of about 30 – 40 hairshafts about one third from the end of the tail. Fold these tightly around the index finger and secure with the thumb.
2. When you have a good grip on the hair, pull quickly.
3. Tape the hair together with masking tape. Masking tape should be used to label the sample with the ear number of the animal. Long ends of hair may be cut off, leaving the roots intact in sample.
4. Referral forms must be completed in full and accompany samples.

Blood samples:

1. Blood samples for determining DNA profiles must be taken in a purple stopper tube, also known as EDTA tube.
2. Approximately 1ml blood is required.
3. Shake the tube thoroughly after taking the sample. Mark the tube clearly.
4. No special transfer requirements is needed.

Samples for the 1/29 translocation:

1. Chromosome studies are needed for identification of the 1/29 translocation.
2. A blood sample in a green stopper tube(Lithium Heparin tube) is required for culture growth.
3. The blood sample for chromosome culturing must not be frozen or transported at very high temperatures.
4. It is important that samples for chromosome studies must be delivered to the laboratory as soon as possible.